## INSTRUCTIONS FOR USE

### General Product Identification Information

*Product supplied in sterile form by ethylene oxide (ETO) according to ISO 11135 - "Preview Sterilization of health-care products -- Ethylene oxide -- Requirements for the development, validation and routine control of a sterilization process for medical devices".*


**Product shelf life: 3 years**

**Technical Name:** Instrument, Manual, Surgical, General Use

**Trade Name:** Meniscus Suture - Out In

**Graphic information that allows the visualization of the final product, description, part number and composition of the material.**

**Table 1.** Commercial models: Graphic Information, Composition, and codes.

CODE	DESCRIPTION	DIMENSION	APPLICATION	RAW-MATERIAL
				
354-100-050	Meniscus Suture - Out in	424 x 30 x 25mm	The Meniscus Suture - Out - In product makes it possible to perform suturing in meniscus operations, both in open and arthroscopic surgeries.	SS ASTM F899, Aluminum NBR ISO 16061, Polyacetal Copolymer NBR ISO16061 and ASTM 1855. Nitinol ASTM 2063

### Identification

The Meniscus Suture - Out In are laser marked with batch number, code, GMReis mark.

### Description of the physical principle and fundamentals of technology, applied for its operation and its action

The Meniscus Suture - Out In is used for meniscus suture surgery and can be used both in open field surgeries and by arthroscopy.

### Support materials that accompany the product for its assembly and protection

All additional information regarding how to use is described in the product usage instruction contained within the product packaging

### How Supplied

The model is supplied sterile using ethylene oxide, duly identified, that is, laser marked with its code, batch number and GMReis logo.

	<h1 style="text-align: center;">INSTRUCTIONS FOR USE</h1> <h2 style="text-align: center;">Meniscus Suture - Out In</h2>		<b>MENOUTIN</b>	
		<b>Sheet:</b>	<b>2 of 5</b>	

Its form of commercial presentation is packaging in double surgical grade paper packaging and packed in a paper box, plasticized, printed, which receives a seal with a Sensitive Indicator of Sterilization by ETO, label, and adhesive security seal, constituting the secondary packaging. Finally, the paper box is wrapped in a polyolefin, heat-shrinkable film, constituting the product's tertiary packaging. The product comes with 5 traceability labels and a manual on how the user can obtain the Instructions for Use of the product, in the non-printed format, at no additional cost, including shipping. This information is compatible with the product labeling. GMReis recommends single use.

AUT. ME: 1.02.477.0  
 CNPJ: 06.045.898/0001-19  
 I.E.: 244.342.283.119  
 Razão Social: G. M. dos Reis Ind. e Com. Ltda  
 Av. Pierre Simon de Laplace, 620, Lote 3 Quadra F  
 CEP: 13045-220 - Campinas/SP - Condomínio Techno Park  
 Telefone: +55 19 3765-8600  
 E-mail: [sac@gmreis.com.br](mailto:sac@gmreis.com.br)

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**Figura 2: Manual that comes with the product.**

It is the responsibility of the medical team and the hospital (the depositories of patient information and records) to attach the label to the surgery description sheet or section of the medical record. GMReis recommends that the patient is the depository of this information, receiving a certificate with the attachment of one of the traceability adhesive labels.

### Special care for storage, conservation, handling and/or transport and associated risks

The product must be stored, handled, and transported in such a way as to prevent any damage or alteration to its characteristics and packaging.

The medical product should be handled with all care, in order to avoid sudden shocks, falls and other risks and / or imperfections that affect the quality of the same and also the safety of the user.

The handling of the product should be done exclusively by medical-hospital professionals, properly trained, qualified and familiar with the technique and procedures involved.

	<b>INSTRUCTIONS FOR USE</b>  Meniscus Suture - Out In		<b>MENOUTIN</b>  Sheet: 3 of 5	
---	---	--	--------------------------------------	--

The effects of vibration, shocks, faulty seating, and improper stacking during transport should be avoided.  
GM REIS RECOMMENDS SINGLE USE.

**If the product is stored, handled and/or transported outside the specific conditions above may pose risks to the procedure and the patient.**

#### **Special precautions, restrictions, warnings, and cautions**

Products that for some reason were not used during the surgical procedure and had their packaging opened, or suffered damage and/or falls cannot be reused and/or reprocessed, and must be discarded in the hospital according to medical-hospital procedure and/or legislation site. If the product presents any abnormality in its operation, a new unit must be opened and put into use.

The person responsible for the procedure must send the non-conforming product to GMReis, duly decontaminated, packaged with batch, product code and with the declaration of the non-compliance that occurred.

#### **Indication, purpose or use for which the product is intended**

##### **Indication**

The Meniscus Suture - Out In is indicated for meniscus suture surgery and can be used both in open field surgeries and by arthroscopy.

##### **Contraindications**

Use of the product is contraindicated for any procedure other than surgical ligament reconstruction procedures.

#### **Characteristics associated with the integrity of the package, way of handling the product and opening the packages in order to ensure sterility**

If any of the packaging is violated or damaged, all components must be discarded, even if they are in perfect condition.

In surgical-grade paper packaging, the integrity of the closure can influence the shelf life of the products packaged in them, since failures in heat sealing constitute entry points for microorganisms, air and/or moisture or even for loss of gaseous mixtures and organic vapors. Therefore, problems with the integrity of the closure can lead to aseptically packaged products being contaminated by the entry of microorganisms.

A closure is considered intact when it ensures the product, at least, the same protection that is conferred by the packaging material.

The integrity of the packaging must be preserved during transport, handling, storage, and distribution of the product.

This product will be supplied sterile using ethylene oxide, therefore, it must be kept in its original packaging until the time of use, in accordance with hospital medical asepsis procedures.

**Restrictions and instructions on packaging damage**

The product must be scrapped and disposed of according to the scrapping and disposal procedure under the following circumstances:

- Violated packaging security seal;
- Physical damage to any of the packages;
- The unreadable label or identification code;
- Shelf life is expired.

**Disuse and disposal**

Before disposal, all components must be unused with the aid of moulders, refans, any other instruments or ways, until it is clearly identified that the product is unsuitable for its use, and there is no possibility of it being reused. After use, the product should be discarded according to the hospital area procedure.

**Additional information about the procedure prior to use of the product**

As the product is supplied sterile, there is no need for any preliminary procedure to use such as cleaning, sterilization and decontamination of the product.

**Physician-specific guidance on reporting adverse events and technical complaints**

If the product presents adverse events not reported in the instruction of use or has technical complaints about the product, the surgeon should immediately contact the manufacturer through the Customer Service (SAC) of GM Reis, in addition to notifying the competent health authority.

To ensure traceability of the product, the surgeon should proceed according to the "product traceability procedure". Product traceability is ensured through the 5 traceability labels provided inside the package, together with the product, as described in the item "Product Traceability Procedure".

**Customer Complaint**

If the medical product presents a specific unpredictable risk, is out of its specifications or is generating any dissatisfaction, notify GM Reis Customer Service (SAC) directly. The product should be forwarded clean and packed in a plastic bag, duly identified and with the description of non-conformity to the following address:

**Av. Pierre Simon de Laplace, 600 - Lote 3 - Quadra F - Block 9677 - ZIP Code 13069-320 - Condominium Technopark - Campinas, SP, Brazil or notify directly at Tel.: (0XX19) 3283-9020; E-mail: [sac@gmreis.com.br](mailto:sac@gmreis.com.br).**

	<b>INSTRUCTIONS FOR USE</b> Meniscus Suture - Out In		<b>MENOUTIN</b> Sheet: 5 of 5	
---	---	--	----------------------------------	--

## Symbology of Labelling

The graphic symbols used in the labeling are in accordance with the NBR ISO 15223 Standard as follows:

SYMBOLS	DESCRIPTIONS	SYMBOLS	DESCRIPTIONS
	Date of manufacture		Single-Use Product "Do not reuse"
	Use-by date		Do not re-sterilize
	Do not use if the package is damaged.		Caution, consult accompanying documents "Attention See Instructions for Use"
	Manufacturer		Batch code
	Sterilized by ethylene oxide		Keep away from light
	Keep dry		Prescription only

## ALERT TO THE CUSTOMER

These Instructions for Use are made available in unprinted format, through the manufacturer's e-mail address: <http://www.gmreis.com.br/produtos/IFU>. Always check the correlation of the Version of the Instructions for Use with the version reported on the label. The instructions for use provided will always be in accordance with the latest version in force. If the user is interested, the Instructions for Use may be made available in printed format at no additional cost. Request for free by e-mail: [sac@gmreis.com.br](mailto:sac@gmreis.com.br).

GM dos Reis Indústria e Comércio Ltda.  
Pierre Simon de Laplace Avenue, nº 600 - Lot 3 - Block F - 9677  
Nova Aparecida Neighborhood  
TECHNOPARK - ZIP Code: 13069-320 – Campinas – SP - BRAZIL  
Operating permit – AFE nº 1.02.477 - 0  
C.N.P.J/M.F: 60.040.599/0001-19 / I.E: 244.342.283.119  
E-mail: [gmreis@gmreis.com.br](mailto:gmreis@gmreis.com.br)  
Phone.: (0XX19) 3765-9900  
Issuing date: 06/24/2022

### Technical and Legal Responsible Qualified:

Geraldo Marins dos Reis Júnior  
CREA - SP N° 0682127536

Revisions History
Rev. 00 06/24/2022 – Notification
Rev. 01 03/01/2023 – Device Listing in USA