


MICRO EASYPASS

SUTURE PASSERS

ENGLISH 



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2022


MICRO EASYPASS

The Micro EasyPass suture passers were developed to facilitate surgical procedures, enabling the transport of sutures through soft tissue or bone tunnels.

The Micro EasyPass consists of a metal cannula and a nitinol wire loop, supplied in a sterile, single-use package, that could be used to pass the Stitch GMReis USP 2-0 to 2 suture.

The variety of models, straight and curved, allows the surgeon to use the Micro EasyPass more suitable for each type of procedure.

Ergonomic handle facilitates the manipulation



High strength nitinol wire loop
Ø0.25 x 560.0 mm

Five options: straight and curved
(tight, small, large, and extra large)



Compatible with Stitch GMReis
UPS 2 - 0 -2 suture



EXTRA LARGE
CURVE

LARGE
CURVE

SMALL
CURVE

TIGHT
CURVE

STRAIGHT

CODE	MODEL
335-600	Micro EasyPass - Micro Suture Passer Extra Large Curve
335-300	Micro EasyPass - Micro Suture Passer Large Curve
335-200	Micro EasyPass - Micro Suture Passer Small Curve
335-99	Micro EasyPass - Micro Suture Passer Tight Curve
335-100	Micro EasyPass - Micro Suture Passer Straight

The Micro EasyPass suture passers are supplied in sterile single-use package, and come with a nitinol wire loop Ø0.25 x 560.0 mm

Stitch - Surgical Sutures and Tapes

High-resistance braided non-absorbable surgical sutures and tapes, made with UHMWPE (ultra high molecular weight polyethylene) polymer, Biocompatible and High resistance against: abrasion, breaking, stretching, loosening or loosening knots.



Micro EasyPass suture passers are compatible with Stitch USP 2-0 and 2, which are sold separately.



Fig.: All surgical suture and suture tapes are supplied with needle.

CODE	MEASURES	COLOR	NEEDLE
333-2-036-AB3817CI	USP 2-0 - 45 cm	Blue and White	3/8 - 17 mm - Taper Point
333-239-AB1226CI	UPS 2 - 90 cm	Blue and White	1/2 - 26 mm - Taper Point
333-239-PB1226CI	USP 2 - 90 cm	Black and White	1/2 - 26 mm - Taper Point



HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone : +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br