

PRESSLOCK

IMPACTED KNOTLESS PEEK ANCHORS

ENGLISH 



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2023

PRESSLOCK ANCHORS Ø2.5 X 8.0 MM

Impacted knotless PEEK Anchors, with closed eyelet, with diameter and length suitable for the treatment of extremities soft tissue injuries: foot, hand and wrist.



Anchor suitable for use with GMReis #2-0 Stitch Suture.



Fig.: TFCC reinsertion with Ø2.5 x 8.0 mm Presslock Anchor and # 2-0 Stitch Suture.



Impaction locking device.

Presslock Anchor provided Cable-mounted insert.


Suture passer device.

PRESSLOCK KNOTLESS IMPACTED PEEK ANCHOR WITH CLOSED EYELET

CODE	Ø	LENGTH
320-25080	2.5 mm	8.0 mm

PRESSLOCK ANCHORS Ø2.9 X 15.5 MM

Impacted knotless PEEK anchors, with closed eyelet, with diameter and length suitable for the treatment of soft tissue injuries in regions such as the shoulder and ankle.



Anchor suitable for use with GMReis # 2 Stitch Suture.





Fig.: Bankart injury treated with Ø2.9 x 15.5 mm Presslock anchors and # 2 Stitch Sutures.



Impaction locking device

Suture passer device.

PRESSLOCK KNOTLESS IMPACTED PEEK ANCHOR WITH CLOSED EYELET


CODE	Ø	LENGTH
320-29155	2.9 mm	15.5 mm

PRESSLOCK HIP ANCHORS Ø2.9 X 15.5 MM

Knotless impacted PEEK anchors with closed eyelet, diameter and length suitable to hip labral injury treatment, with long inserter and specially developed instruments for arthroscopic surgery.




Fig.: Hip labral injury treated with GMReis Ø2.9 x 15.5 mm Presslock Hip and # 2 Stitch Sutures




PRESSLOCK HIP KNOTLESS IMPACTED PEEK ANCHOR WITH CLOSED EYELET

CODE	Ø	LENGTH
320-29155-SFPK-HIP	2.9 mm	10.5 / 15.5 mm



Ø2.5 x 8.0 mm PRESSLOCK INSTRUMENTAL

CODE	DESCRIPTION
320-CA-24	Presslock Cannula Ø2.4 mm
320-TR-24	Presslock Trocar Tip Ø2.4 mm
320-110-25	Presslock Drill Ø2.4 mm



Ø2.9 x 15.5 mm PRESSLOCK INSTRUMENTAL

CODE	DESCRIPTION
320-CA-29	Presslock Cannula Ø2.9 mm
320-TR-29	Presslock Trocar Tip Ø2.9 mm
320-110-29	Presslock Drill Ø2.9 mm




Ø2.9 X 15.5 MM PRESSLOCK HIP INSTRUMENTAL

CODE	DESCRIPTION
320-CA-29-HIP	Presslock Hip Cannula Ø2.9 mm
320-TR-29-HIP	Presslock Hip Trocar Tip Ø2.9 mm
320-110-29-HIP	Presslock Hip Drill Ø2.9 mm

STITCH - SURGICAL SUTURE AND TAPES

Presslock anchors do not come with surgical suture and tapes; for Combined use, Stitch Sutures and Tapes are recommended.



High strength sutures and tapes no absorbables, manufactured with UHMWPE polymer (high molecular weight polyethylene), biocompatible and with high resistance against: abrasion, breakage, stretching, loosening or release of the Knots.


Pic.: All Stitch Suture and Tapes are supplied needled.

UHMWPE WIRE WITH NEEDLE

CODE	MEASURES	COLOR	NEEDLE
333-2-036-AB3817CI	USP 2-0 – 45 cm	Blue and White	3/8 – 17 mm - Cylindrical Tip
333-239-AB1226CI	USP 2 – 90 cm	Blue and White	1/2 – 26 mm - Cylindrical Tip
333-239-PB1226CI	USP 2 – 90 cm	Black and White	1/2 – 26 mm - Cylindrical Tip

UHMWPE TAPE WITH NEEDLE

CODE	MEASURES	COLOR	NEEDLE
333-20-AB	2.0 mm x 100 cm	Blue and White	1/2 - 26 mm - Cylindrical Tip
333-20-B	2.0 mm x 100 cm	Black and White	1/2 - 26 mm - Cylindrical Tip





HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br