


# ARTROM

ARTHROEREISIS SCREW

ENGLISH 



**GMREIS**

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة


2024

# ARTHROEREISIS SCREW - ARTROM

Artrom is a conical and cannulated subtalar implant made of titanium for talocalcaneal arthroereisis; developed for treatment of various disorders related to flat foot, which result in subtalar joint instability and excessive pronation.

The conical shape provides a perfect anatomical fit of Artrom in relation to sinus tarsi region; threaded body reduces irritation of the surrounding tissues and prevents implant migration.

*T30 hexalobular connection provides better coupling and torque transfer between implant and instruments, ensuring safety in placement and removal of Artrom.*




## Indications

- Severe hindfoot pronation;
- Calcaneus valgus;
- Exaggerated plantar flexion of talus;
- Painful flat foot;
- As an adjunct to treatment of tarsal coalitions;
- Instability of subtalar joint;
- Flat foot acquired from adult due to posterior tibial dysfunction and,
- Paralytic flat foot.

## ARTROM ARTHROEREISIS SCREW







CODE	Ø	LENGHT	COLOR
241-70-12	7.0 mm	12.0 mm	
241-80-14*	8.0 mm	14.0 mm	
241-90-14	9.0 mm	14.0 mm	
241-100-14	10.0 mm	14.0 mm	
241-110-16	11.0 mm	16.0 mm	
241-120-16	12.0 mm	16.0 mm	

\*Consult availability, sale upon prior request.



CODE	DESCRIPTION
241-32-EST	Percutaneous Guide for Arthroereisis



CODE	DESCRIPTION	QUANTITY	COLOR
241-200	Screwdriver - Artrom	01	-
241-300	Artrom Probe	01	-
241-7012-T	Arthroereisis Template Ø7 mm	01	
241-8014-T*	Arthroereisis Template Ø8 mm	01	
241-9014-T	Arthroereisis Template Ø9 mm	01	
241-10014-T	Arthroereisis Template Ø10 mm	01	
241-11016-T	Arthroereisis Template Ø11 mm	01	
241-12016-T	Arthroereisis Template Ø12 mm	01	
241-1000	Artrom Case	01	

\*Consult availability, sale upon prior request.



**HEADQUARTER | GMREIS**

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park  
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone : +55 (19) 3765 9900  
marketingdigital@gmreis.com.br

 GMReisBrasil

[www.gmreis.com.br](http://www.gmreis.com.br)