


MICRO EASYPASS

SUTURE PASSERS

ENGLISH 



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

2024

MICRO EASYPASS

The Micro EasyPass suture passers were developed to facilitate surgical procedures, enabling the transport of sutures through soft tissue or bone tunnels.

The Micro EasyPass consists of a metal cannula and a nitinol wire loop, supplied in a sterile, single-use package, that could be used to pass the Stitch GMReis USP 2-0 to 2 suture.

The variety of models, straight and curved, allows the surgeon to use the Micro EasyPass more suitable for each type of procedure.



High strength nitinol wire loop Ø0.25 x 560.0 mm.

Five options: straight and curved (tight, small, large, and extra large)



Compatible with Stitch GMReis UPS 2 - 0 -2 suture

Ergonomic handle facilitates the manipulation




| CODE | MODEL |
|---------|--|
| 335-600 | Micro EasyPass - Micro Suture Passer Extra Large Curve |
| 335-300 | Micro EasyPass - Micro Suture Passer Large Curve |
| 335-200 | Micro EasyPass - Micro Suture Passer Small Curve |
| 335-99 | Micro EasyPass - Micro Suture Passer Tight Curve |
| 335-100 | Micro EasyPass - Micro Suture Passer Straight |

The Micro EasyPass suture passers are supplied in sterile single-use package, and come with a nitinol wire loop Ø0.25 x 560.0 mm

STITCH - SURGICAL SUTURES AND TAPES

High-resistance braided non-absorbable surgical sutures and tapes, made with UHMWPE (ultra high molecular weight polyethylene) polymer, Biocompatible and High resistance against: abrasion, breaking, stretching, loosening or loosening knots.



Micro EasyPass suture passers are compatible with Stitch USP 2-0 and 2, which are sold separately.



Fig.: All surgical suture and suture tapes are supplied with needle.

| CODE | MEASURES | COLOR | NEEDLE |
|--------------------|-----------------|-----------------|---------------------------|
| 333-2-036-AB3817CI | USP 2-0 - 45 cm | Blue and White | 3/8 - 17 mm - Taper Point |
| 333-239-AB1226CI | UPS 2 - 90 cm | Blue and White | 1/2 - 26 mm - Taper Point |
| 333-239-PB1226CI | USP 2 - 90 cm | Black and White | 1/2 - 26 mm - Taper Point |



GMREIS

Qualidade para Vida

Calidad para Vida Quality for Life الجودة للحياة

HEADQUARTER | GMREIS

Pierre Simon de Laplace, 600 - Lote 3 Quadra F - Techno Park
ZIP CODE: 13069-320 | Campinas / SP - Brazil | Phone .: +55 (19) 3765 9900
marketingdigital@gmreis.com.br

 GMReisBrasil

www.gmreis.com.br